



CNK1 Monoclonal Antibody

Catalog No	YP-mAb-03778
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CNKSR1
Immunogen	The antiserum was produced against synthesized peptide derived from human CNKR1. AA range:491-540
Specificity	CNK1 Monoclonal Antibody detects endogenous levels of CNK1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNK;CNK homolog protein 1;CNK1;CNKSR1;Connector enhancer of KSR 1; Connector enhancer of KSR like;hCNK1;KSR;Connector enhancer of kinase suppressor of ras 1;Connector enhancer of KSR-like
Observed Band	80kD
Calculated Molecular Weight	80kD
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein .
Tissue Specificity	Brain,Epithelium,Kidney,Liver,Neuroepithelium,Place
Function	function:May function as an adapter protein or regulator of Ras signaling pathways.,PTM:Phosphorylated on tyrosine.,similarity:Belongs to the CNKSR family.,similarity:Contains 1 CRIC domain.,similarity:Contains 1 PDZ (DHR) domain.,similarity:Contains 1 PH domain.,similarity:Contains 1 SAM (sterile alpha motif) domain.,subunit:Interacts with RHO and RALGDS.,
Background	This gene encodes a protein containing several motifs involved in protein-protein interaction, including PDZ, PH (Pleckstrin homology), and SAM (sterile alpha motif) domains. The encoded protein acts as a scaffold component for receptor tyrosine kinase signaling and may mediate crosstalk between different signaling



pathways. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CNK1 Monoclonal Antibody