




DACA-1 Monoclonal Antibody

Catalog No	YP-mAb-03811
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	YTHDF1
Immunogen	The antiserum was produced against synthesized peptide derived from human YTHDF1. AA range:1-50
Specificity	DACA-1 Monoclonal Antibody detects endogenous levels of DACA-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf21;DACA-1;DF1;YTH domain-containing family protein 1;YTHDF1;YTH domain family protein 1;Dermatomyositis associated with cancer putative autoantigen 1
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	Cytoplasm . Cytoplasm, P-body . Cytoplasm, Stress granule .
Tissue Specificity	Cervix carcinoma,Colon,Ileal mucosa,Lymph,Melanoma,Placenta
Function	similarity:Contains 1 YTH domain.,
Background	similarity:Contains 1 YTH domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DACA-1 Monoclonal Antibody