



eIF4E3 Monoclonal Antibody

Catalog No	YP-mAb-03849
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EIF4E3
Immunogen	The antiserum was produced against synthesized peptide derived from human EIF4E3. AA range:141-190
Specificity	eIF4E3 Monoclonal Antibody detects endogenous levels of eIF4E3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	eIF-4E3;eIF-4E type 3;eIF4E type 3;eIF4E 3;eIF 4E3;eIF4E3;Eukaryotic translation initiation factor 4E type 3;eIF4E-3
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	cytoplasm,cytosol,mRNA cap binding complex,
Tissue Specificity	Brain,
Function	function:Recognizes and binds the 7-methylguanosine-containing mRNA cap during an early step in the initiation of protein synthesis. May act as an inhibitor of EIF4E1 activity.,similarity:Belongs to the eukaryotic initiation factor 4E family., subunit:eIF4F is a multi-subunit complex, the composition of which varies with external and internal environmental conditions. It is composed of at least eIF4A, eIF4E and eIF4G (By similarity). EIF4E3 interacts with EIF4G1, but not with EIF4EBP1, EIF4EBP2 and EIF4EBP3.,
Background	EIF4E3 belongs to the EIF4E family of translational initiation factors that interact with the 5-prime cap structure of mRNA and recruit mRNA to the ribosome (Joshi et al., 2004 [PubMed 15153109]).[supplied by OMIM, Mar 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using eIF4E3 Monoclonal Antibody