



Eps8L2 Monoclonal Antibody

Catalog No	YP-mAb-03857
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EPS8L2
Immunogen	The antiserum was produced against synthesized peptide derived from human ES8L2. AA range:1-50
Specificity	Eps8L2 Monoclonal Antibody detects endogenous levels of Eps8L2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EPS8 like 2;EPS8 like protein 2;EPS8 related protein 2;EPS8L2;EPS8R2; Epidermal growth factor receptor kinase substrate 8-like protein 2;PP13181; EPS8-like protein 2;Epidermal growth factor receptor pathway substrate 8-related protein 2;EPS8-related protein 2
Observed Band	81kD
Calculated Molecular Weight	81kD
Cell Pathway	Cytoplasm . Cell projection, stereocilium . Localizes at the tips of the stereocilia of the inner and outer hair cells. .
Tissue Specificity	Detected in fibroblasts and placenta.
Function	function:Stimulates guanine exchange activity of SOS1. May play a role in membrane ruffling and remodeling of the actin cytoskeleton.,similarity:Belongs to the EPS8 family.,similarity:Contains 1 PID domain.,similarity:Contains 1 SH3 domain.,subunit:Interacts with ABI1. Part of a complex that contains SOS1, ABI1 and EPS8L2. Associates with F-actin.,tissue specificity:Detected in fibroblasts and placenta.,



Background

This gene encodes a member of the EPS8 gene family. The encoded protein, like other members of the family, is thought to link growth factor stimulation to actin organization, generating functional redundancy in the pathways that regulate actin cytoskeletal remodeling. [provided by RefSeq, Dec 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Eps8L2 Monoclonal Antibody