



HRF Monoclonal Antibody

Catalog No	YP-mAb-03913
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPT1
Immunogen	The antiserum was produced against synthesized peptide derived from human TCTP. AA range:71-120
Specificity	HRF Monoclonal Antibody detects endogenous levels of HRF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fortilin;Histamine-releasing factor;HRF;p23;TCTP;TPT1;Translationally-controlled tumor protein
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Found in several healthy and tumoral cells including erythrocytes, hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells, gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.
Function	function:Involved in calcium binding and microtubule stabilization., similarity:Belongs to the TCTP family.,subunit:Interacts with STEAP3.,tissue specificity:Found in several healthy and tumoral cells including erythrocytes, hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells, gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.,



Background

function:Involved in calcium binding and microtubule stabilization.,
similarity:Belongs to the TCTP family.,subunit:Interacts with STEAP3.,tissue
specificity:Found in several healthy and tumoral cells including erythrocytes,
hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells,
gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected
in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HRF Monoclonal Antibody