



HSP70L Monoclonal Antibody

Catalog No	YP-mAb-03915
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSPA1L/HSPA1A
Immunogen	The antiserum was produced against synthesized peptide derived from human HSPA1L. AA range:221-270
Specificity	HSP70 Monoclonal Antibody detects endogenous levels of HSP70 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Heat shock 7 kDa protein 14;Heat shock 7kDa protein 14;Heat shock protein HSP6;Heat shock protein hsp7 related protein;HSP6;HSP7 4;HSP7 like protein 1;HSP7-like protein 1;HSP7L1;HSP7E_HUMAN;HSPA 14;HSPA14;MGC13199;HSP70L;HSPA1L/HSPA1A;Heat shock 70 kDa protein 1-like/Heat shock 70 kDa protein 1A/1B;HSPA1L;Heat shock 70 kDa protein 1-like;Heat shock 70 kDa protein 1L;Heat shock 70 kDa protein 1-Hom;HSP70-Hom;HSPA1A;HSPA1;HSPA1B;Heat shock 70 kDa protein 1A/1B;Heat shock 70 kDa protein 1/2;HSP70-1/HSP70-2;HSP70.1/HSP70.2
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	zona pellucida receptor complex,nucleoplasm,mitochondrial matrix,cytosol,COP9 signalosome,cell body,blood microparticle,
Tissue Specificity	Expressed in spermatids.
Function	function:In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within organelles. These chaperones participate in all these processes through their ability to recognize nonnative



conformations of other proteins. They bind extended peptide segments with a net hydrophobic character exposed by polypeptides during translation and membrane translocation, or following stress-induced damage.,induction:Not induced by heat shock.,similarity:Belongs to the heat shock protein 70 family.,tissue specificity:Expressed in spermatids.,

Background

This gene encodes a 70kDa heat shock protein. In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the folding of newly translated proteins in the cytosol and in organelles. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which also encode isoforms of the 70kDa heat shock protein. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HSP70L Monoclonal Antibody