



HSP90A Monoclonal Antibody

Catalog No	YP-mAb-03925
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSP90AA1
Immunogen	The antiserum was produced against synthesized peptide derived from human HSP90A. AA range:683-732
Specificity	HSP90A Monoclonal Antibody detects endogenous levels of HSP90A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSP90AA1;EC:3.6.4.10;Heat shock protein family C member 1;Heat shock protein HSP 90 alpha;Heat shock protein HSP 90-alpha;HSP 90;HSP90A;HSPC1;HSPCA;Heat shock 86 kDa;HSP 86;HSP86;Renal carcinoma antigen NY-REN-38
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus . Cytoplasm . Melanosome . Cell membrane . Mitochondrion . Identified by mass spectrometry in melanosome fractions from stage I to stage IV.
Tissue Specificity	Adult brain,Brain,Brain cortex,Breast,Cajal-Retzius
Function	caution:Despite classification as a pseudogene, the existence of this protein is supported by unambiguous mass spectrometry evidence.,function:Molecular chaperone.,function:Molecular chaperone. Has ATPase activity.,similarity:Belongs to the heat shock protein 90 family.,subcellular location:Identified by mass spectrometry in melanosome fractions from stage I to stage IV.,subunit:Homodimer. Interacts with AHSA1, SMYD3 and TOM34. Interacts with FNIP1 and HSF1.,



Background

The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HSP90A Monoclonal Antibody