




Mob3C Monoclonal Antibody

Catalog No	YP-mAb-03967
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOB3C
Immunogen	The antiserum was produced against synthesized peptide derived from human MOL2C. AA range:81-130
Specificity	Mob3C Monoclonal Antibody detects endogenous levels of Mob3C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MOB kinase activator 3C;Mob1 homolog 2C;MOB1E;MOB3C;Mps One Binder kinase activator like 2C;MOBKL2C;Mps one binder kinase activator-like 2C
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	
Tissue Specificity	Lung,Spleen,
Function	function:May regulate the activity of kinases.,similarity:Belongs to the MOB1/phocein family.,
Background	The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Mob3C Monoclonal Antibody