



MOV10L1 Monoclonal Antibody

Catalog No	YP-mAb-03968
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MOV10L1
Immunogen	The antiserum was produced against synthesized peptide derived from human MOV10L1. AA range:318-367
Specificity	MOV10L1 Monoclonal Antibody detects endogenous levels of MOV10L1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cardiac helicase activated by MEF2 protein;Cardiac specific RNA helicase; CHAMP;Csm;M1L1_HUMAN;Moloney leukemia virus 1 like protein 1;Moloney leukemia virus 1-like protein 1;MOV1 like 1;MOV1 like protein 1;MOV1-like protein 1;MOV1L1;Moloney leukemia virus 1 like 1;homolog (mouse);Putative helicase Mov1l1;MOV10L1;Putative helicase Mov10l1;Moloney leukemia virus 10-like protein 1;MOV10-like protein 1
Observed Band	140kD
Calculated Molecular Weight	140kD
Cell Pathway	Cytoplasm . Component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. .
Tissue Specificity	[Isoform 1]: Specifically expressed in testis.
Function	alternative products:Additional isoforms seem to exist,function:Putative RNA helicase. Isoform 1 may play a role in male germ cell development., similarity:Belongs to the DNA2/NAM7 helicase family. SDE3 subfamily.,tissue specificity:Isoform 1 is specifically expressed in testis.,



Background

Mov10 RISC complex RNA helicase like 1(MOV10L1) Homo sapiens This gene is similar to a mouse gene that encodes a putative RNA helicase and shows testis-specific expression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MOV10L1 Monoclonal Antibody