



MRP-L10 Monoclonal Antibody

Catalog No	YP-mAb-03973
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL10
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL10. AA range:211-260
Specificity	MRP-L10 Monoclonal Antibody detects endogenous levels of MRP-L10 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRP-L10;MRP L8;MRP L10;Large ribosomal subunit protein uL10m;L8mt;MRPL10;39S ribosomal protein L10 mitochondrial;MRPL8;RPML8;39S ribosomal protein L10;mitochondrial;L10mt;39S ribosomal protein L8;MRP-L8;mitochondrial
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Pancreas,Skin,
Function	similarity:Belongs to the ribosomal protein L10P family.,subunit:Component of the mitochondrial ribosome large subunit (39S) which comprises a 16S rRNA and about 50 distinct proteins.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian



mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Sequence analysis identified three transcript variants that encode two different isoforms. A pseudogene corresponding to this gene is found on chromosome 5q. [provided by RefSeq, Nov 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-L10 Monoclonal Antibody