




MRP-L14 Monoclonal Antibody

Catalog No	YP-mAb-03977
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL14
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL14. AA range:96-145
Specificity	MRP-L14 Monoclonal Antibody detects endogenous levels of MRP-L14 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRP-L14;MRP L32;MRP L14;Large ribosomal subunit protein uL14m;L32mt;MRPL14;39S ribosomal protein L14 mitochondrial;MRPL32;RPML32;39S ribosomal protein L14;mitochondrial;L14mt;39S ribosomal protein L32;MRP-L32;mitochondrial
Observed Band	16kD
Calculated Molecular Weight	16kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Pancreas,
Function	similarity:Belongs to the ribosomal protein L14P family.,
Background	This nuclear gene encodes a protein component of the 39S subunit of the mitochondrial ribosome. A pseudogene of this gene is found on chromosome 17. [provided by RefSeq, Jan 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using MRP-L14 Monoclonal Antibody