



MRP-L18 Monoclonal Antibody


Catalog No	YP-mAb-03981
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL18
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL18. AA range:1-50
Specificity	MRP-L18 Monoclonal Antibody detects endogenous levels of MRP-L18 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRP-L18;MRPL 18;MRP L18;Large ribosomal subunit protein uL18m;L18mt;MRPL18;39S ribosomal protein L18 mitochondrial;HSPC071;39S ribosomal protein L18;mitochondrial
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Brain, Eye, Umbilical cord blood,
Function	similarity:Belongs to the ribosomal protein L18P family.,
Background	This nuclear gene encodes a protein component of the larger 39S subunit of mitochondrial ribosome. This protein may also aid in the import of nuclear-encoded 5S rRNA into mitochondria. Alternative splicing results in multiple transcript variants, most of which are not predicted to encode a protein. A pseudogene of this gene is found on chromosome 16. [provided by RefSeq, Jan 2016],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using MRP-L18 Monoclonal Antibody