



MRP-L20 Monoclonal Antibody

Catalog No	YP-mAb-03982
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL20
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL20. AA range:100-149
Specificity	MRP-L20 Monoclonal Antibody detects endogenous levels of MRP-L20 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L20mt;Large ribosomal subunit protein bL20m;MRP L20;MRP-L20;MRPL20;39S ribosomal protein L20 mitochondrial;39S ribosomal protein L20;mitochondrial
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Bone marrow,Lymph,Mammary gland,Whole embryo,
Function	similarity:Belongs to the ribosomal protein L20P family.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. A



pseudogene corresponding to this gene is found on chromosome 21q. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2016],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-L20 Monoclonal Antibody