



MRP-L40 Monoclonal Antibody

Catalog No	YP-mAb-03992
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL40
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL40. AA range:101-150
Specificity	MRP-L40 Monoclonal Antibody detects endogenous levels of MRP-L40 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NLVCF;MRP-L40;MRP L40;MRP L22;Large ribosomal subunit protein mL40;MRPL40;39S ribosomal protein L40 mitochondrial;URIM;39S ribosomal protein L40;mitochondrial;L40mt;Nuclear localization signal-containing protein deleted in velocardiofacial syndrome;Up-regulated in metastasis
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Ubiquitous.
Function	tissue specificity:Ubiquitous.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly



in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Deletions in this gene may contribute to the etiology of velo-cardio-facial syndrome and DiGeorge syndrome. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-L40 Monoclonal Antibody