



MRP-S21 Monoclonal Antibody

Catalog No	YP-mAb-04010
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPS21
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPS21. AA range:38-87
Specificity	MRP-S21 Monoclonal Antibody detects endogenous levels of MRP-S21 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small ribosomal subunit protein bS21m;S21mt;RPMS21;MRP-S21;MRP S21;MRPS21;28S ribosomal protein S21 mitochondrial;MDS016;28S ribosomal protein S21;mitochondrial
Observed Band	11kD
Calculated Molecular Weight	11kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Brain,Cerebellum,Hematopoietic stem cell,Lung,Testis,Uterus,
Function	similarity:Belongs to the ribosomal protein S21P family.,subunit:Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA.



Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that belongs to the ribosomal protein S21P family. Pseudogenes corresponding to this gene are found on chromosomes 1p, 1q, 9p, 10p, 10q, 16q, and 17q. Available sequence data analyses identified splice varian


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-S21 Monoclonal Antibody