



MRP-S34 Monoclonal Antibody

Catalog No	YP-mAb-04018
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPS34
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPS34. AA range:169-218
Specificity	MRP-S34 Monoclonal Antibody detects endogenous levels of MRP-S34 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small ribosomal subunit protein mS34;S34mt;MRP-S34;MRPS12;MRP S34;MRPS34;28S ribosomal protein S34 mitochondrial;28S ribosomal protein S34;mitochondrial
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Eye,
Function	subunit:Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly



in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-S34 Monoclonal Antibody