



MRP-S36 Monoclonal Antibody

Catalog No	YP-mAb-04020
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPS36
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPS36. AA range:4-53
Specificity	MRP-S36 Monoclonal Antibody detects endogenous levels of MRP-S36 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S36mt;MRP S36;KGD4;DC47;Alpha-ketoglutarate dehydrogenase subunit 4; MRP-S36;MRPS36;28S ribosomal protein S36 mitochondrial;28S ribosomal protein S36;mitochondrial
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Dendritic cell,Kidney,Prostate,
Function	subunit:Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. The mitochondrial ribosome (mitoribosome) consists of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly



in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. Pseudogenes corresponding to this gene are found on chromosomes 3p, 4q, 8p, 11q, 12q, and 20p. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRP-S36 Monoclonal Antibody