



NEDL2 Monoclonal Antibody

Catalog No	YP-mAb-04035
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HECW2
Immunogen	The antiserum was produced against synthesized peptide derived from human HECW2. AA range:481-530
Specificity	NEDL2 Monoclonal Antibody detects endogenous levels of NEDL2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin-protein ligase HECW2;EC:2.3.2.26;HECT;C2 and WW domain-containing protein 2;NEDL2;HECW2;KIAA1301;NEDD4-like E3 ubiquitin-protein ligase 2
Observed Band	140-150kD
Calculated Molecular Weight	140-150kD
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton, spindle .
Tissue Specificity	Predominantly expressed in adult brain, lung and heart.
Function	function:E3 ubiquitin-protein ligase that mediates ubiquitination of TP73. Acts to stabilize TP73 and enhance activation of transcription by TP73.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 C2 domain.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,similarity:Contains 2 WW domains.,subunit:Interacts with TP73.,tissue specificity:Predominantly expressed in adult brain, lung and heart.,
Background	function:E3 ubiquitin-protein ligase that mediates ubiquitination of TP73. Acts to stabilize TP73 and enhance activation of transcription by TP73.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 C2 domain.,



similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,
similarity:Contains 2 WW domains.,subunit:Interacts with TP73.,tissue
specificity:Predominantly expressed in adult brain, lung and heart.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using NEDL2 Monoclonal Antibody