



Peroxin 14 Monoclonal Antibody

Catalog No	YP-mAb-04069
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PEX14
Immunogen	The antiserum was produced against synthesized peptide derived from human PEX14. AA range:117-166
Specificity	Peroxin 14 Monoclonal Antibody detects endogenous levels of Peroxin 14 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAPP2;Peroxin 14;Peroxin-14;PEX;Pex14p;PEX14;Peroxisomal membrane protein PEX14;PTS1 receptor-docking protein;Peroxisomal membrane anchor protein PEX14
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Peroxisome membrane ; Peripheral membrane protein ; Cytoplasmic side .
Tissue Specificity	Brain,Cerebellum,Epithelium,Muscle,Placenta,Testis,
Function	disease:Defects in PEX14 are a cause of Zellweger syndrome (ZWS) [MIM:214100]. ZWS is a fatal peroxisome biogenesis disorder characterized by dysmorphic facial features, hepatomegaly, ocular abnormalities, renal cysts, hearing impairment, profound psychomotor retardation, severe hypotonia and neonatal seizures. Death occurs within the first year of life.,disease:Defects in PEX14 are the cause of peroxisome biogenesis disorder complementation group K (PBD-CGK) [MIM:601791]. PBD refers to a group of peroxisomal disorders arising from a failure of protein import into the peroxisomal membrane or matrix. The PBD group is comprised of four disorders: Zellweger syndrome



(ZWS), neonatal adrenoleukodystrophy (NALD), infantile Refsum disease (IRD), and classical rhizomelic chondrodysplasia punctata (RCDP). ZWS, NALD and IRD are distinct from RCDP and constitute a clinical continuum of overlapping

Background

peroxisomal biogenesis factor 14(PEX14) Homo sapiens This gene encodes an essential component of the peroxisomal import machinery. The protein is integrated into peroxisome membranes with its C-terminus exposed to the cytosol, and interacts with the cytosolic receptor for proteins containing a PTS1 peroxisomal targeting signal. The protein also functions as a transcriptional corepressor and interacts with a histone deacetylase. A mutation in this gene results in one form of Zellweger syndrome. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Peroxin 14 Monoclonal Antibody