



PGBD1 Monoclonal Antibody

Catalog No	YP-mAb-04075
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PGBD1
Immunogen	The antiserum was produced against synthesized peptide derived from human PGBD1. AA range:241-290
Specificity	PGBD1 Monoclonal Antibody detects endogenous levels of PGBD1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SCAND4;PiggyBac transposable element-derived protein 1;PGBD 1;hucep-4; HUCEP4;PGBD1;Cerebral protein 4
Observed Band	92kD
Calculated Molecular Weight	92kD
Cell Pathway	nucleus,membrane,
Tissue Specificity	Brain,Liver,
Function	similarity:Contains 1 SCAN box domain.,
Background	The piggyBac family of proteins, found in diverse animals, are transposases related to the transposase of the canonical piggyBac transposon from the moth, Trichoplusia ni. This family also includes genes in several genomes, including human, that appear to have been derived from the piggyBac transposons. This gene belongs to the subfamily of piggyBac transposable element derived (PGBD) genes. The PGBD proteins appear to be novel, with no obvious relationship to other transposases, or other known protein families. This gene product is specifically expressed in the brain, however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding the



same protein.[provided by RefSeq, May 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PGBD1 Monoclonal Antibody