




Pmp24 Monoclonal Antibody

Catalog No	YP-mAb-04095
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PXMP4
Immunogen	The antiserum was produced against synthesized peptide derived from human PXMP4. AA range:15-64
Specificity	Pmp24 Monoclonal Antibody detects endogenous levels of Pmp24 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Pmp24;PXMP4;Peroxisomal membrane protein 4;24 kDa peroxisomal intrinsic membrane protein
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Peroxisome membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in normal prostate epithelial cells, and androgen-sensitive prostate adenocarcinoma cells. Not expressed in androgen-insensitive prostate adenocarcinoma cells.
Function	similarity:Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit:Interacts with PEX19.,
Background	similarity:Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit:Interacts with PEX19.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Pmp24 Monoclonal Antibody