




PWWP2B Monoclonal Antibody

Catalog No	YP-mAb-04106
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PWWP2B
Immunogen	The antiserum was produced against synthesized peptide derived from human PWWP2B. AA range:268-317
Specificity	PWWP2B Monoclonal Antibody detects endogenous levels of PWWP2B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PWWP Domain Containing 2;PWWP Domain Containing 2B;PWWP domain-containing protein 2B;PWWP2B;PWWP2
Observed Band	64kD
Calculated Molecular Weight	64kD
Cell Pathway	nucleoplasm,
Tissue Specificity	Placenta,Uterus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 PWWP domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 PWWP domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PWWP2B Monoclonal Antibody