




RFPL4A Monoclonal Antibody

Catalog No	YP-mAb-04117
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RFPL4A
Immunogen	The antiserum was produced against synthesized peptide derived from human RFPL4A. AA range:155-204
Specificity	RFPL4A Monoclonal Antibody detects endogenous levels of RFPL4A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ret finger protein like 4;Ret finger protein-like 4A;RFPL4;RFPL4A;RNF210;RING finger protein 210
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	
Background	The ret finger protein-like 4 gene, Rfpl4, encodes a putative E3 ubiquitin-protein ligase expressed in adult germ cells.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RFPL4A Monoclonal Antibody