




Ribosomal Protein L26L Monoclonal Antibody

Catalog No	YP-mAb-04135
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL26L1
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL26L. AA range:96-145
Specificity	Ribosomal Protein L26L Monoclonal Antibody detects endogenous levels of Ribosomal Protein L26L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	60S ribosomal protein L26-like 1;Large ribosomal subunit protein uL24-like 1; Ribosomal Protein L26L;RPL26L1;RPL26P1
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	cytosol,large ribosomal subunit,cytosolic large ribosomal subunit,extracellular exosome,
Tissue Specificity	Lung,Umbilical cord blood,
Function	similarity:Belongs to the ribosomal protein L24P family.,
Background	This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Dec 2015],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Ribosomal Protein L26L Monoclonal Antibody