




Ribosomal Protein S27L Monoclonal Antibody

Catalog No	YP-mAb-04166
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPS27L
Immunogen	The antiserum was produced against synthesized peptide derived from human RPS27L. AA range:1-50
Specificity	Ribosomal Protein S27L Monoclonal Antibody detects endogenous levels of Ribosomal Protein S27L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small ribosomal subunit protein eS27-like;RPS;Ribosomal protein eS27-like; Ribosomal Protein S27L;RPS27L;40S ribosomal protein S27-like
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	nucleus,nucleoplasm,ribosome,cytosolic small ribosomal subunit,
Tissue Specificity	Kidney,Umbilical cord blood,
Function	cofactor: Binds 1 zinc ion per subunit .,similarity:Belongs to the ribosomal protein S27e family.,
Background	This gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27, which suggests the encoded protein may be a component of the 40S ribosomal subunit. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Ribosomal Protein S27L Monoclonal Antibody