




SAA4 Monoclonal Antibody

Catalog No	YP-mAb-04184
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human SAA4 (NP_006503.2).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C SAA;CSAA;SAA4;Serum amyloid A 4 protein;serum amyloid A4;Serum amyloid A-4 protein;Constitutively expressed serum amyloid A protein;C-SAA
Observed Band	15kDa
Calculated Molecular Weight	14kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Predicted to be involved in acute-phase response. Located in extracellular exosome.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using SAA4 Mouse mAb (A16428) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 90s.