



SCCA2 Monoclonal Antibody

Catalog No	YP-mAb-04186
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERPINB4
Immunogen	The antiserum was produced against synthesized peptide derived from human SERPINB4. AA range:261-310
Specificity	SCCA2 Monoclonal Antibody detects endogenous levels of SCCA2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LEUPIN;Peptidase inhibitor 11;PI 11;PI11;SCCA;SCCA 2;SCCA1;SCCA2;Serpine B4;SERPINB3;SERPINB4;PI-11;Squamous cell carcinoma antigen 2;SCCA-2
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Cytoplasm . Seems to also be secreted in plasma by cancerous cells but at a low level.
Tissue Specificity	Squamous cells.
Function	function:May act as a protease inhibitor to modulate the host immune response against tumor cells.,similarity:Belongs to the serpin family.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,subcellular location:Seems to also be secreted in plasma by cancerous cells but at a low level.,tissue specificity:Squamous cells.,
Background	function:May act as a protease inhibitor to modulate the host immune response against tumor cells.,similarity:Belongs to the serpin family.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,subcellular location:Seems to also be



secreted in plasma by cancerous cells but at a low level.,tissue specificity:Squamous cells.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SCCA2 Monoclonal Antibody