



TIS11D Monoclonal Antibody

Catalog No	YP-mAb-04246
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZFP36L2
Immunogen	The antiserum was produced against synthesized peptide derived from human TISD. AA range:207-256
Specificity	TIS11D Monoclonal Antibody detects endogenous levels of TIS11D protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRF2;Butyrate response factor 2;EGF response factor 2;EGF-response factor 2;ERF 2;TIS11D;ZFP36L2;Zinc finger protein 36 C3H1 type-like 2;ERF2;RNF162C;Zinc finger protein 36;C3H1 type-like 2;ZFP36-like 2;ERF-2;Protein TIS11D
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Nucleus . Cytoplasm . Shuttles between the nucleus and the cytoplasm in a XPO1/CRM1-dependent manner. .
Tissue Specificity	Expressed mainly in the basal epidermal layer, weakly in the suprabasal epidermal layers (PubMed:27182009). Expressed in epidermal keratinocytes (at protein level) (PubMed:27182009).
Function	function:Probable regulatory protein involved in regulating the response to growth factors. RNA-binding protein that binds to 5'UUUUUUUUU-3' core sequence. Binds to the class II AU-rich element (ARE) in the 3'-UTR of target mRNAs and promotes their deadenylation and degradation.,similarity:Contains 2 C3H1-type zinc fingers.,



Background

This gene is a member of the TIS11 family of early response genes. Family members are induced by various agonists such as the phorbol ester TPA and the polypeptide mitogen EGF. The encoded protein contains a distinguishing putative zinc finger domain with a repeating cys-his motif. This putative nuclear transcription factor most likely functions in regulating the response to growth factors. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TIS11D Monoclonal Antibody