




TRIM59 Monoclonal Antibody

Catalog No	YP-mAb-04256
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM59
Immunogen	The antiserum was produced against synthesized peptide derived from human TRIM59. AA range:191-240
Specificity	TRIM59 Monoclonal Antibody detects endogenous levels of TRIM59 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRF1;RING finger protein 104;RNF104;TRIM57;TRIM59;tripartite motif containing 59;TSBF1;Tumor;Tripartite motif-containing protein 59;Tumor suppressor TSBF-1
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	Human testis,
Function	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,
Background	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TRIM59 Monoclonal Antibody