




ZIP9 Monoclonal Antibody

Catalog No	YP-mAb-04276
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC39A9
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC39A9. AA range:101-150
Specificity	ZIP9 Monoclonal Antibody detects endogenous levels of ZIP9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC39A9;Zinc transporter ZIP9;ZIP 9;ZIP9;Zrt and Irt like protein 9;Solute carrier family 39 member 9;Zrt- and Irt-like protein 9;ZIP-9
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Coronary artery,Duodenum,Placenta,Skin,
Function	function:May act as a zinc-influx transporter.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.,
Background	function:May act as a zinc-influx transporter.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZIP9 Monoclonal Antibody