



Kallistatin Monoclonal Antibody

Catalog No	YP-mAb-04317
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERPINA4
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human SERPINA4. AA range:161-210
Specificity	Kallistatin Monoclonal Antibody detects endogenous levels of Kallistatin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kallistatin;Serpina4;KAL;Kallikrein inhibitor;KLST;SERPINA4;KST;PI4;Peptidase inhibitor 4;PI-4
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Secreted.
Tissue Specificity	Expressed by the liver and secreted in plasma.
Function	function:Inhibits human amidolytic and kininogenase activities of tissue kallikrein. Inhibition is achieved by formation of an equimolar, heat- and SDS-stable complex between the inhibitor and the enzyme, and generation of a small C-terminal fragment of the inhibitor due to cleavage at the reactive site by tissue kallikrein.,miscellaneous:Heparin blocks kallistatin's complex formation with tissue kallikrein and abolishes its inhibitory effect on tissue kallikrein's activity., PTM:The N-terminus is blocked.,similarity:Belongs to the serpin family., subunit:Monomer and some homodimers.,tissue specificity:Expressed by the liver and secreted in plasma.,



Background

function: Inhibits human amidolytic and kininogenase activities of tissue kallikrein. Inhibition is achieved by formation of an equimolar, heat- and SDS-stable complex between the inhibitor and the enzyme, and generation of a small C-terminal fragment of the inhibitor due to cleavage at the reactive site by tissue kallikrein., miscellaneous: Heparin blocks kallistatin's complex formation with tissue kallikrein and abolishes its inhibitory effect on tissue kallikrein's activity., PTM: The N-terminus is blocked., similarity: Belongs to the serpin family., subunit: Monomer and some homodimers., tissue specificity: Expressed by the liver and secreted in plasma.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Kallistatin Monoclonal Antibody