



Lipocalin-1 Monoclonal Antibody

Catalog No	YP-mAb-04319
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LCN1
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human LCN1. AA range:11-60
Specificity	Lipocalin-1 Monoclonal Antibody detects endogenous levels of Lipocalin-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LCN1;Lipocalin 1;lipocalin 1 (tear prealbumin);PMFA;Tear lipocalin;Tear prealbumin;Tlc;TP;VEG;Lipocalin-1;VEGP;Von Ebner gland protein;VEG protein
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	Secreted.
Tissue Specificity	Mainly expressed in lachrymal and salivary glands. Also expressed in the prostate.
Function	caution:Could be the product of a pseudogene.,function:Could play a role in taste reception. Could be necessary for the concentration and delivery of sapid molecules in the gustatory system. Can bind various ligands, with chemical structures ranging from lipids and retinoids to the macrocyclic antibiotic rifampicin and even to microbial siderophores. Exhibits an extremely wide ligand pocket.,function:May bind a variety of ligands including lipids.,similarity:Belongs to the calycin superfamily. Lipocalin family.,subunit:Homodimer (By similarity). Binds to LMBR1L which may mediate its endocytosis.,tissue specificity:Mainly expressed in lachrymal and salivary glands. Also expressed in the prostate.,



Background

This gene encodes a member of the lipocalin family of small secretory proteins. Lipocalins are extracellular transport proteins that bind to a variety of hydrophobic ligands. The encoded protein is the primary lipid binding protein in tears and is overproduced in response to multiple stimuli including infection and stress. The encoded protein may be a marker for chromosome aneuploidy as well as an autoantigen in Sjogren's syndrome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and two pseudogenes of this gene are also located on the long arm of chromosome 9. [provided by RefSeq, Nov 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Lipocalin-1 Monoclonal Antibody