



Protein C Monoclonal Antibody

Catalog No	YP-mAb-04320
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PROC
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human PROC. AA range:181-230
Specificity	Protein C Monoclonal Antibody detects endogenous levels of Protein C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A1 activator;Cerebroside sulfate activator;Co beta glucosidase;Component C; CSAct;Protein C;PROC;Vitamin K-dependent protein C;Anticoagulant protein C; Autoprothrombin IIA;Blood coagulation factor XIV
Observed Band	52kD
Calculated Molecular Weight	52kD
Cell Pathway	Secreted . Golgi apparatus . Endoplasmic reticulum .
Tissue Specificity	Plasma; synthesized in the liver.
Function	catalytic activity:Degradation of blood coagulation factors Va and VIIIa., disease:Defects in PROC are the cause of protein C deficiency autosomal dominant (ADPROCD) [MIM:176860]. ADPROCD is a cause of hereditary thrombophilia, a hemostatic disorder characterized by impaired regulation of blood coagulation and a tendency to recurrent venous thrombosis. However, many adults with heterozygous disease may be asymptomatic. Individuals with decreased amounts of protein C are classically referred to as having type I protein C deficiency and those with normal amounts of a functionally defective protein as having type II deficiency.,disease:Defects in PROC are the cause of protein C deficiency autosomal recessive (ARPROCD) [MIM:612304].



ARPROCD results in a thrombotic condition that can manifest as a severe neonatal disorder or as a milder disorder with late-onset thrombophilia. The severe form I

Background

This gene encodes a vitamin K-dependent plasma glycoprotein. The encoded protein is cleaved to its activated form by the thrombin-thrombomodulin complex. This activated form contains a serine protease domain and functions in degradation of the activated forms of coagulation factors V and VIII. Mutations in this gene have been associated with thrombophilia due to protein C deficiency, neonatal purpura fulminans, and recurrent venous thrombosis.[provided by RefSeq, Dec 2009],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Protein C Monoclonal Antibody