



DCI Monoclonal Antibody

Catalog No	YP-mAb-04321
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ECI1
Immunogen	Synthesized peptide derived from DCI . at AA range: c-terminal
Specificity	DCI Monoclonal Antibody detects endogenous levels of DCI protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3;2-trans-enoyl-CoA isomerase;D3;D2-enoyl-CoA isomerase;Delta(3);Delta(2)-enoyl-CoA;DCI;ECI1;Enoyl-CoA delta isomerase 1 mitochondrial;Enoyl-CoA delta isomerase 1;mitochondrial;Delta(2)-enoyl-CoA isomerase;Dodecenoyl-CoA isomerase; mitochondrial
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	Expressed in liver (at protein level).
Function	catalytic activity:(3Z)-dodec-3-enoyl-CoA = (2E)-dodec-2-enoyl-CoA., function:Able to isomerize both 3-cis and 3-trans double bonds into the 2-trans form in a range of enoyl-CoA species.,pathway:Lipid metabolism; fatty acid beta-oxidation.,similarity:Belongs to the enoyl-CoA hydratase/isomerase family., subunit:Homotrimer.,
Background	This gene encodes a member of the hydratase/isomerase superfamily. The protein encoded is a key mitochondrial enzyme involved in beta-oxidation of unsaturated fatty acids. It catalyzes the transformation of 3-cis and 3-trans-enoyl-CoA esters arising during the stepwise degradation of cis-, mono-, and



polyunsaturated fatty acids to the 2-trans-enoyl-CoA intermediates. Alternatively spliced transcript variants have been described. [provided by RefSeq, May 2010]


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DCI Monoclonal Antibody