



Mob3A/B Monoclonal Antibody


Catalog No	YP-mAb-04345
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MMOBKL2A/MMOBKL2B
Immunogen	Synthesized peptide derived from Mob3A/B . at AA range: 1-80
Specificity	Mob3A/B Monoclonal Antibody detects endogenous levels of Mob3A/B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MOB LAK;Mob1 homolog 2A;MOB3A;MOBKL2A;moblak;Protein Mob3A;Mob3A/B;MMOBKL2A/MMOBKL2B;Mps one binder kinase activator-like 2A/B;MOB kinase activator 3A;MOB-LAK;Mps one binder kinase activator-like 2A;MOB3B;MOBKL2B;MOB kinase activator 3B;Mob1 homolog 2b;Mps one binder kinase activator-like 2B
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	intracellular,
Tissue Specificity	Cervix,Skin,
Function	function:May regulate the activity of kinases.,sequence caution:Chimera. There is an inclusion of DNA from another gene in position 141.,similarity:Belongs to the MOB1/phocein family.,
Background	function:May regulate the activity of kinases.,sequence caution:Chimera. There is an inclusion of DNA from another gene in position 141.,similarity:Belongs to the MOB1/phocein family.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Mob3A/B Monoclonal Antibody