




FAM3D Monoclonal Antibody

Catalog No	YP-mAb-04349
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM3D
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human FAM3D. AA range:121-170
Specificity	FAM3D Monoclonal Antibody detects endogenous levels of FAM3D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ567/PRO1130;Protein FAM3D;OIT1;EF7;FAM3D
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Secreted .
Tissue Specificity	Abundantly expressed in placenta and weakly expressed in small intestine.
Function	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
Background	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FAM3D Monoclonal Antibody