




TP8L1 Monoclonal Antibody

Catalog No	YP-mAb-05009
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TNFAIP8L1
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	TP8L1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Oxidative stress-regulated gene-beta;Oxy-beta;TIPE1;TNF alpha-induced protein 8-like protein 1;TP8L1;TNFAIP8L1;Tumor necrosis factor alpha-induced protein 8-like protein 1 (TNF alpha-induced protein 8-like protein 1) (TNFAIP8-like protein 1)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Cytoplasm .
Tissue Specificity	High expression detected in most carcinoma cell lines, especially in cells transformed with virus genomes.
Function	similarity:Belongs to the TNFAIP8 family.,
Background	similarity:Belongs to the TNFAIP8 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TP8L1 Monoclonal Antibody