




BARH2 Monoclonal Antibody

Catalog No	YP-mAb-05027
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BARHL2
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	BARH2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BarH like 2 homeobox protein;BarH like homeobox 2;BARHL2;BARH2;BarH-like 2 homeobox protein
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Potential regulator of neural basic helix-loop-helix genes., similarity:Belongs to the BAR homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
Background	function:Potential regulator of neural basic helix-loop-helix genes., similarity:Belongs to the BAR homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using BARH2 Monoclonal Antibody