




F215A Monoclonal Antibody

Catalog No	YP-mAb-05049
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM215A APR2 C17orf88
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	F215A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F215A;FAM215A APR2 C17orf88;Putative protein FAM215A (Putative apoptosis-related protein 2) (APR-2)
Observed Band	12kD
Calculated Molecular Weight	12kD
Cell Pathway	
Tissue Specificity	
Function	apoptosis, cell death, programmed cell death, death,
Background	FAM215A (Family With Sequence Similarity 215 Member A (Non-Protein Coding)) is an RNA Gene, and is affiliated with the non-coding RNA class. It is a Mouse Monoclonal antibody to 6D11 to C17orf88, APR-2, FAM215A, family with sequence similarity 215, member A (non-protein coding), LINC00530. This antibody recognizes Human antigen.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using F215A Monoclonal Antibody