




FGFP3 Monoclonal Antibody

Catalog No	YP-mAb-05074
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FGFBP3 C10orf13 PSEC0101
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	FGFP3 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FGFP3;FGFBP3 C10orf13 PSEC0101;Fibroblast growth factor-binding protein 3 (FGF-BP3) (FGF-binding protein 3) (FGFBP-3)
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,
Function	similarity:Belongs to the fibroblast growth factor-binding protein family., subunit:Interacts with FGF1.,
Background	similarity:Belongs to the fibroblast growth factor-binding protein family., subunit:Interacts with FGF1.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FGFP3 Monoclonal Antibody