




IFN16 Monoclonal Antibody

Catalog No	YP-mAb-05091
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IFNA16
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	IFN16 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IFN16;IFNA16;Interferon alpha-16 (IFN-alpha-16) (Interferon alpha-WA)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Secreted.
Tissue Specificity	
Function	function:Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.,similarity:Belongs to the alpha/beta interferon family.,
Background	function:Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.,similarity:Belongs to the alpha/beta interferon family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using IFN16 Monoclonal Antibody