




CAB39 Monoclonal Antibody

Catalog No	YP-mAb-05399
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CAB39 MO25 CGI-66
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CAB39 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein Mo25;MO25;CGI 66;calcium binding protein 39;CAB39;CAB39 MO25 CGI-66;Calcium-binding protein 39 (MO25alpha) (Protein Mo25)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Duodenum,Hypothalamus,Testis,
Function	function:Together with the STE20-related adaptor-alpha (STRAD alpha) pseudo kinase, forms a regulatory complex caMABLE of stimulating the activity of STK11.,similarity:Belongs to the Mo25 family.,
Background	function:Together with the STE20-related adaptor-alpha (STRAD alpha) pseudo kinase, forms a regulatory complex caMABLE of stimulating the activity of STK11.,similarity:Belongs to the Mo25 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CAB39 Monoclonal Antibody