




ZC3H1 Monoclonal Antibody

Catalog No	YP-mAb-05484
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZFC3H1 CCDC131 KIAA0546 PSRC2
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZC3H1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC131;CSRC2;KIAA0546;PSRC2;ZC3H1;ZFC3H1 CCDC131 KIAA0546 PSRC2;Zinc finger C3H1 domain-containing protein (Coiled-coil domain-containing protein 131) (Proline/serine-rich coiled-coil protein 2)
Observed Band	218kD
Calculated Molecular Weight	218kD
Cell Pathway	Nucleus . Excluded from the nucleolus. .
Tissue Specificity	Blood,Brain,Epithelium,Lymph,Placenta,Thalamus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 5 TPR repeats., similarity:Contains 6 HAT repeats.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 5 TPR repeats., similarity:Contains 6 HAT repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZC3H1 Monoclonal Antibody