




LDH6A Monoclonal Antibody

Catalog No	YP-mAb-05697
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LDHAL6A LDHL2
Immunogen	Synthesized peptide derived from human protein . at AA range: 50-130
Specificity	LDH6A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L-lactate dehydrogenase A-like 6A;LDHL2;LDHA-like protein 6A;EC:1.1.1.27; LDH6A;LDHAL6A LDHL2;L-lactate dehydrogenase A-like 6A (EC 1.1.1.27)
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis-specific.
Function	catalytic activity:(S)-lactate + NAD(+) = pyruvate + NADH., pathway:Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.,similarity:Belongs to the LDH/MDH superfamily., similarity:Belongs to the LDH/MDH superfamily. LDH family.,
Background	catalytic activity:(S)-lactate + NAD(+) = pyruvate + NADH., pathway:Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.,similarity:Belongs to the LDH/MDH superfamily., similarity:Belongs to the LDH/MDH superfamily. LDH family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using LDH6A Monoclonal Antibody