



RAB15 Monoclonal Antibody


Catalog No	YP-mAb-06037
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RAB15
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	RAB15 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Rab 15;RAB15;RAB15 member RAS onocogene family;member RAS onocogene family;RAB15_HUMAN;Ras related protein Rab 15;Ras related protein Rab15;Ras-related protein Rab-15
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	Brain,Fetal brain,Placenta,
Function	function:May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:The GTP bound form of RAB15 interacts with REP15.,
Background	function:May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:The GTP bound form of RAB15 interacts with REP15.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RAB15 Monoclonal Antibody