




K0087 Monoclonal Antibody

Catalog No	YP-mAb-06373
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIAA0087 HA1002
Immunogen	Synthesized peptide derived from human protein . at AA range: 20-100
Specificity	K0087 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0087;K0087;KIAA0087 HA1002;Uncharacterized protein KIAA0087
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	
Tissue Specificity	Bone marrow,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using K0087 Monoclonal Antibody