




SFT2A Monoclonal Antibody

Catalog No	YP-mAb-06392
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SFT2D1 C6orf83
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SFT2A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SFT2A;SFT2D1 C6orf83;Vesicle transport protein SFT2A (SFT2 domain-containing protein 1) (pRGR1)
Observed Band	17kD
Calculated Molecular Weight	17kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Choriocarcinoma,
Function	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SFT2A Monoclonal Antibody