




HSPB9 Monoclonal Antibody

Catalog No	YP-mAb-06620
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSPB9
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	HSPB9 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSPB9;Heat shock protein beta-9 (HspB9) (Cancer/testis antigen 51) (CT51)
Observed Band	17kD
Calculated Molecular Weight	17kD
Cell Pathway	Cytoplasm . Nucleus . Translocates to nuclear foci during heat shock.
Tissue Specificity	Testis specific.
Function	similarity:Belongs to the small heat shock protein (HSP20) family.,tissue specificity:Testis specific.,
Background	similarity:Belongs to the small heat shock protein (HSP20) family.,tissue specificity:Testis specific.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using HSPB9 Monoclonal Antibody