




CREG2 Monoclonal Antibody

Catalog No	YP-mAb-07202
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CREG2
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	CREG2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein CREG2;CREG2
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Secreted .
Tissue Specificity	Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.
Function	PTM:It is not sure whether N-glycosylation is on Asn-165 and/or Asn-166., similarity:Belongs to the CREG family.,tissue specificity:Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.,
Background	PTM:It is not sure whether N-glycosylation is on Asn-165 and/or Asn-166., similarity:Belongs to the CREG family.,tissue specificity:Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CREG2 Monoclonal Antibody