




COLGALT1 mouse mAb

Catalog No	YP-mAb-07890
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLT25D1 PSEC0241
Immunogen	Synthesized peptide derived from human COLGALT1 AA range: 194-244
Specificity	This antibody detects endogenous levels of COLGALT1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.04% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PSEC0241;EC:2.4.1.50;Collagen beta(1-O)galactosyltransferase 1;COLGALT1; ColGalT 1;GLT25D1 PSEC0241;Procollagen galactosyltransferase 1 (EC 2.4.1.50) (Glycosyltransferase 25 family member 1) (Hydroxylysine galactosyltransferase 1)
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Endoplasmic reticulum lumen . Colocalized with PLOD3 and mannose binding lectin/MBL2. .
Tissue Specificity	Ubiquitous with higher levels in placenta, heart, lung and spleen.
Function	similarity:Belongs to the glycosyltransferase 25 family.,
Background	The protein encoded by this gene is one of two enzymes that transfers galactose moieties to hydroxylysine residues of collagen and mannose binding lectin. This gene is constitutively expressed and encodes a soluble protein that localizes to the endoplasmic reticulum. [provided by RefSeq, Dec 2015],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using COLGALT1 mouse mAb